

Exporting Sequence of Selected Annotations

Open the *Sequence View* with a document that contains annotations. A good candidate here could be any file in Genbank format with both sequence and annotations. Select a single or several annotations or annotation groups in the *Annotation editor*, click the right mouse button to open the context menu and select the *Export* *Export sequence of selected annotations* item:

The screenshot displays a sequence viewer interface with a context menu open over a selected region. The menu items include: Find qualifier..., Invert annotation selection, Go to position... (Ctrl+G), Select sequence region... (Ctrl+A), New annotation... (Ctrl+N), Copy/Paste, Select, Add, Analyze, Align, Cloning, Fetch sequences from remote database, Export (highlighted), Edit, Remove, Rulers..., and Disable 'CDS' highlighting. The 'Export' menu item has opened a sub-menu with three options: Save circular view as image..., Export selected sequence region..., and Export sequence of selected annotations... (highlighted). The background shows sequence tracks with annotations like 'CDS' and 'codon start'.

The *Export Sequence of Selected Annotations* dialog will appear which is similar to the *Export Selected Sequences* dialog described [here](#).