


Creating New Project

A *project* store links to the data files, cross-file data associations, and visualization settings.

Below is a description of how to create a new project manually. Note that if you have no project created when opening a file with a sequence, an alignment or any other biological data, a new anonymous project is created automatically.

To create a new project select the *File* *New project* menu or click the *New project* button on the main toolbar.

The dialog will appear:



The dialog box is titled "Create New Project". It features two input fields: "Project name" containing "New Project" and "Save project to file" containing "/home/..|project.uprj". The "Save project to file" field includes a file selection icon (three dots). At the bottom, there are three buttons: "Help", "Cancel", and "Create".

Here you need to specify the visual name for the project and the directory and file to store it.

After you click the *Create* button the [Project View](#) window is opened.