

# Exporting Sequence Image

Use a sequence toolbar *Export image* button to save a screenshot of the sequence:



The Export Image dialog will appear where you should set the name, location, export settings, and format of the picture:

### Export Image

**Export settings**

**Area**

- Currently viewed
- Zoomed annotations
- Sequence details

**Export to file**

File name:  ...

Format:  ▼

UGENE supports export to the BMP, JPG, PNG, TIFF and SVG image formats. You can export currently viewed, zoomed annotations or sequence details areas. Also, you can export the whole sequence or custom region. Use the *Region* settings to do it.