

Options Panel

The *Options Panel* is available in the *Sequence View*, *Alignment Editor* and in the *Assembly Browser*. By default, it is closed.

To open a tab of the Options Panel click on the corresponding icon at the right side of a *Sequence View*, *Alignment Editor*, *Sanger Editor* or *Assembly Browser* window. To close the tab click again on the tab icon.

More detailed information about different *Options Panel* tabs can be found in the following chapters:

- *Options Panel in Sequence View*
 - [Search in Sequence](#)
 - [Highlighting Annotations](#)
 - [Information about Sequence](#) (Statistics)
 - [In Silico PCR](#)
 - [Circular View Settings](#)
- *Options Panel in Alignment Editor*
 - [Consensus](#)
 - [Searching in Alignment](#)
 - [Highlighting Alignment](#)
 - [Pairwise Alignment](#)
 - [Tree Settings](#)
 - [Export Consensus](#)
 - [Statistics](#)
- *Options Panel in Sanger Editor*
 - [Alignment Statistics](#)
 - [Export Chromatogram Consensus](#)
 - [Sanger Reads Settings](#)
- *Options Panel in Assembly Browser*
 - [Navigation in Assembly Browser](#)
 - [Assembly Browser Settings](#)
 - [Assembly Statistic](#)