Remote BLASTing

The workflow sample, described below, allows one to do remote queries to the NCBI BLAST database to search for homologous nucleotide sequences for multiple input sequences at the same time.

As the result of the BLAST each input sequence is annotated with the "blast result" annotations. These annotations are used to fetch the corresponding homologous sequences from the NCBI database based on the identifiers specified in the "blast result" annotations. The output homologous sequences and the original sequences, annotated by BLAST, are grouped by folders.



If you haven't used the workflow samples in UGENE before, look at the "How to Use Sample Workflows" section of the documentation.

Workflow Sample Location

The workflow sample "Remote BLASTing" can be found in the "Scenarios" section of the Workflow Designer samples.

Workflow Image

The opened workflow looks as follows:



Workflow Wizard

The wizard has 3 pages.

1. Input Sequence(s) Page: On this page you must input at least one nucleotide sequence.



- my_sequence2.fa

 \oslash

2. <u>Remote Nucleotide BLAST Page:</u> Here you can optionally modify parameters that should be used for the remote BLAST queries. For example, you can select the search database, correct the e-value and set the maximum number of results (i.e. "Max hits"). The "Megablast" option, applied by default, specifies to optimize the search for high similar sequences only. Selecting it decreases the search time, but some less similar results could be skipped by the search in this case. Note that the "Megablast" option is also applied by default in the NCBI BLAST web interface.

Remote BLASTing Wizard			? ×
	Remote Nucleotide BLAST		
	Remote BLAST Parame Database		
	Expected value	10.000000	×
	Max hits Megablast	10	
		True	
	Additional Show additional parameters		+
Defaults		< Back Next >	Cancel

There are also some additional parameters. Description of them can be found in the Remote BLAST workflow element chapter of the documentation.

Additional Hide additional paramete	rs -
Gap costs	02
Match scores	1-3
Short sequence	False
Entrez query	
Annotate as	

3. <u>Output Files Page:</u> this is an informational page. It states that this workflow has predefined names of the output files. For each input sequence the workflow outputs:

- "orig_with_blast.gb" file: the file contains the input sequence itself and the "blast result" annotations;
 "homologous.gb" file: the file contains the found homologous sequences loaded from the NCBI by identifiers, specified in the BLAST results.

The results on the hard drive are grouped by folders (see below).

The wizard page looks as follows:



Workflow Result

The workflow output files are shown in the dashboard as follows:

UGENE - [Workflow Designer - Remote BLASTing]					- • ×	
					_ 8 ×	
				\$ <mark>0</mark>	To Workflow Designer	
/ Remote BLASTing 1 💥						
Output Files		Workflow Task				
File	Producer	Time: 00:01:14				
orig_with_blast.gb 🗸	Write Original Sequence(s) with BLAS					
homologous.gb 🗸	Write Homologous Sequences	The workflow tas	k has been finished	successfully!		
orig_with_blast.gb 🗸	Write Original Sequence(s) with BLAS	Common Statistics				
homologous.gb 🗸	Write Homologous Sequences	Element	Elapsed time	Output message	æ	
		Remote BLAST	00:01:09	2		
		Fetch Sequences by ID from	00:00:01	20		
		Read Sequence(s)	00:00:00	2		
		Write Original Sequence(s) w	00:00:00	0		
		Write Homologous Sequences	00:00:00	0		
2: Tasks 🛄 <u>3</u> : Log					No active tasks 💡 🚺	

Each file can be opened in the UGENE Sequence View by clicking on the corresponding link in the dashboard.

On the hard drive the output is grouped by folders with the names of the input sequences. For example, for the input sequences specified above, the output hierarchy will be the following:

- my_sequence1.fa folder with files:
- my_sequence1.ta tolder with files:
 orig_with_blast.gb
 homologous.gb
 my_sequence2.ta folder with files:
 orig_with_blast.gb
 homologous.gb