

Highlighting Annotations

To configure settings of annotation names go to the *Annotation Highlighting* tab in the *Options Panel*.

By default the tab shows annotations names of the opened *Sequence View*.

NC_00136

15002k2.5k3k3.5k4k4.5k5k5833

5' terminusteature (2)3' terminust

15002k2.5k3k3.5k4k4.5k5k5833

3875←[1125 bp]→4999

DHPLD★HGAFNAM

TIL★TDMATHSTPC

GPPSSRLTWRIQRH

GACCATCCTCTAGACTGACATGGCGCATTCAACGCCCATG

385638603865387038753880388538903894

CTGGTAGGAGATCTGACTGTACCGCGTAAGTTGCGGTAC

GDELSVH★RM★RW

V★MR★V★SMACEVGH

SWGRSQCPANLALM

NameType

▶CDS

▶CDS

▶CDS

▶CDS

▶comment (0,1)

▶misc_feature (0,2)

▶misc_feature

▶misc_feature

▶source (0,1)

CDS

CDS

CDS

CDS

Misc

Misc

Annotations Highlighting

Select an annotation name:

| Annotation | Color |
|--------------|-------|
| BadI | |
| BaeGI | |
| BaeI | |
| CDS | |
| SacII | |
| comment | |
| misc_feature | |
| source | |

Show all annotation names

Configure the annotations:

☒ Show annotations

☐ Show on translation

☒ Show value of qualifier:

label,note

←

→

Help

If you want to see all annotation names, click the *Show all annotation names* link. The *Previous annotation* and *Next annotation* buttons seek to the previous or to the next annotation of the view correspondingly.

Find below information about annotations names' properties that you can configure.

- [Annotations Color](#)
- [Annotations Visability](#)
- [Show on Translation](#)
- [Captions on Annotations](#)