

# Selecting Sequence Region

When you [are selecting a region of a sequence](#) e.g in the *Sequence zoom view* the corresponding region on the 3D structure is being highlighted while the rest regions of the 3D structure are being shaded.

To configure the color of a region selected open the *Settings* dialog (press the *Settings* item in the 3D Structure Viewer context menu or in the *Display* menu on the toolbar to do it), press the *Set selection color* button and select a color in the dialog appeared.

To adjust the shading drag the *Unselected regions shading* slider in the *Settings* dialog.

