


Exporting Sequence Image

Use a sequence toolbar *Export image* button to save a screenshot of the sequence:



The Export Image dialog will appear where you should set name, location, export settings and format of the picture:

A screenshot of the 'Export Image' dialog box. The dialog has a title bar with a close button. It is divided into three main sections. The first section, 'Export settings', contains two sub-sections: 'Area' with three radio buttons ('Currently viewed', 'Zoomed annotations', and 'Sequence details', with 'Sequence details' selected) and 'Region' with a dropdown menu set to 'Selected region' and two input boxes containing '31' and '75' separated by a hyphen. The second section, 'Export to file', contains a 'File name' field with a browse button ('...'), a 'Format' dropdown set to 'png', and a 'Quality' slider set to '80'. The third section contains three buttons: 'Help', 'Cancel', and 'Export'.

UGENE supports export to the BMP, JPEG, JPG, PNG, PPM, TIF, TIFF, XBM, XPM and SVG image formats. You can export currently viewed, zoomed annotations or sequence detailes areas. Also you can export whole sequence or custom region. Use the *Region* settings to do it.