

Quality Filter

This workflow filters sequences with quality \geq than parameter "quality" and writes result in file in FASTQ format.



How to Use This Sample


If you haven't used the workflow samples in UGENE before, look at the ["How to Use Sample Workflows"](#) section of the documentation.

Workflow Sample Location

The workflow sample "Quality Filter" can be found in the "Custom Elements" section of the Workflow Designer samples.

Workflow Image


The workflow looks as follows:




Workflow Wizard

The wizard has 2 pages.

1. Input sequence(s): On this page you must input sequence(s).



2. Quality Filter: On this page you can modify quality filter and output settings.



The following parameters are available:

| | |
|-----------------------|--|
| Minimum quality value | Minimum quality value |
| Result FASTQ file | Location of output data file. If this attribute is set, slot "Location" in port will not be used. |
| Accumulate results | Accumulate all incoming data in one file or create separate files for each input. In the latter case, an incremental numerical suffix is added to the file name. |