

Find Correct Primer Pairs Element

Find correct primer pairs, which consist of valid primers without dimers.

Parameters in GUI

| Parameter | Description | Default value |
|--------------------|---------------------------------|---------------|
| Output report file | Path to the report output file. | |

Parameters in Workflow File

Type: find-primers

| Parameter | Parameter in the GUI | Type |
|-------------|----------------------|---------------|
| output-file | Output report file | <i>string</i> |

Input/Output Ports

The element has 1 *input port*.

Name in GUI: Input sequences

Name in Workflow File: in-sequence

Slots:

| Slot In GUI | Slot in Workflow File | Type |
|-------------|-----------------------|----------|
| Sequence | sequence | sequence |