

Aligning with ClustalW

Task Name: align-clustalw

Multiple sequence alignment with ClustalW.

 ClustalW is used as an *external tool* and must be installed on your system.

Parameters:

toolpath — path to the ClustalW executable. By default, the path specified in the *Application Settings* is applied. [String, Optional, Default: "default"]

tmpdir — directory for temporary files. [String, Optional]

in — semicolon-separated list of input files. [String, Required]

out — output file. [String, Required]

format — format of the output file. [String, Optional]

Example:

```
ugene align-clustalw --in=C0I.aln --out=C0I.sto --format=stockholm
```