

Find Correct Primer Pairs Element

Find correct primer pairs, which consist of valid primers without dimers.

Parameters in GUI

Parameter	Description	Default value
Output report file	Path to the report output file.	

Parameters in Workflow File

Type: find-primers

Parameter	Parameter in the GUI	Type
output-file	Output report file	string

Input/Output Ports

The element has 1 *input port*:

Name in GUI: Input sequences

Name in Workflow File: in-sequence

Slots:

Slot In GUI	Slot in Workflow File	Type
Sequence	sequence	<i>sequence</i>