


Export Chromatogram Consensus

To export consensus sequence use the *Consensus* tab of the *Options Panel*:



The screenshot shows the 'Consensus' tab of an 'Options Panel'. It features a 'Consensus mode' section with a 'Consensus type' dropdown set to 'Strict' and a 'Threshold' slider set to 100%. Below this is an 'Export consensus' section with a file path input 'Aligned reads_consensus_3.txt', a 'File format' dropdown set to 'Plain text', a 'Keep gaps' checkbox, and an 'Export' button. The panel has a light beige background with a red-shaded area for the export options. A vertical toolbar on the right contains icons for a gear, a document, and a list.

The following parameters are available:

Export to file - here you need to select path for the output file.

File format - format for the output file.

Keep gaps - set this parameter if the result consensus must keep gaps.

When you click on the *Export* button the consensus sequence will be exported into selected output file.