

Export Chromatogram Consensus

To export consensus sequence use the *Consensus* tab of the *Options Panel*:

Consensus

▼ Consensus mode

Consensus type:
Strict

Threshold:
100%

▼ Export consensus

Export to file
Aligned reads_consensus_3.txt ...

File format
Plain text

☐ Keep gaps

Export

The following parameters are available:

Export to file - here you need to select path for the output file.

File format - format for the output file.

Keep gaps - set this parameter if the result consensus must keep gaps.

When you click on the *Export* button the consensus sequence will be exported into selected output file.