

Highlighting Annotations

To configure settings of annotation types go to the *Annotation Highlighting* tab in the *Options Panel*.

By default the tab shows annotations types of the opened *Sequence View*.

NC_001363 [dna]

15001k1.5k2k2.5k3k3.5k4k4.5k5k5833

CDS (4)
NP_040335.1NP_5977NP_5977NP_040336.1NP
5' terminal
source
3' terminal

15001k1.5k2k2.5k3k3.5k4k4.5k5k5833

F I * S P E K G G N E R P H P * V A
L F S L Q K K G G M K D P T R R W Q
Y L V S R K R G E * K T P P V G G
TTATTTAGTCTCCAGAAAAAGGGGGGAATGAAGACCCACCCGTAGGTGGCA
52195225523052355240524552505255526052655271
AATAAATCAGAGGTCTTTTCCCCCTTACTTTCTGGGGTGGGCATCCACCGT
* K T E L F L P S H F V G G T P P L
K N L R W F F P P I F S G V R L H C
I * D G S F P P F S L G W G Y T A

Name	Value
Auto-annotations [murine.gb NC_001363]	
NC_001363 features [murine.gb]	
CDS (0, 4)	
CDS	1042..2658
CDS	join(2970..3413,3412..3873)
CDS	3875..4999
CDS	5048..5203
misc_feature (0, 2)	
misc_feature	2..590
misc_feature	5245..5833

Annotations Highlighting

Select an annotation type:

Annotation	Color
CDS	
misc_feature	
source	

Show all annotation types

Configure the annotation type:

☒ Show annotations of this type

☐ Show on translation

☒ Show value of qualifier:

←

→

Help

If you want to see all annotation types, click the *Show all annotation types* link. The *Previous annotation* and *Next annotation* buttons seek to the previous or to the next annotation of the view correspondingly.

Find below information about annotations types' properties that you can configure.

- [Annotations Color](#)
- [Annotations Visability](#)
- [Show on Translation](#)
- [Captions on Annotations](#)