Exporting Chromatogram Data

Open, for example, the \$UGENE/data/samples/SCF/90-JRI-07.srf file. In the Project View context menu there is Export chromatogram to SCF item:



After clicking on the item, the Export chromatogram file dialog will appear:

atogram File
C:/work/ugene/data/samples/SCF/90-JRI-07_copy.scf
SCF 🔹
Complemented
to the project
Export Cancel Help

Check the Reversed and Complemented options if you want to create a reverse and complement chromatogram. Press the Export button.

The exported file will be opened in the Sequence View.