Get Amino Translations of a Sequence

The workflow takes a nucleotide sequence as input and returns two files: translations of the sequence and translations of the complement sequence.

How to Use This Sample


If you haven't used the workflow samples in UGENE before, look at the "How to Use Sample Workflows" section of the documentation.

Workflow Sample Location

The workflow sample "Get Amino Translations of a Sequence" can be found in the "Scwnarios" section of the Workflow Designer samples.

Workflow Image

The opened workflow looks as follows:



Workflow Wizard

The wizard has 2 pages.

1. Input sequence(s): On this page you must input sequence(s).

Get Amino Translations Wizard	8 ×
Input sequence(s)	
Defaults Next >	Cancel Help

2. <u>Output Data:</u> On this page you can modify result files settings.

Get Amino Translations Wizard				8 ×
	Output data Result files Direct translations Compl. translations	direct_translations.fa		
Defaults	< Bac	k Apply	Cancel Run	Help