

Exporting Alignment as Image

To export an alignment as an image click the *Export as image* button on the editor toolbar or call the *Export->Export as image* context menu item.



The Export Image dialog will appear where you should set name, location, export settings and format of the picture:

Export Image

Export settings

Region Whole alignment Select

☐ Include sequences' names

☐ Include consensus

☒ Include ruler

Export to file

File name ...

Format PNG

Help Cancel Export

UGENE supports export to the BMP, CUR, DDS, ICNS, JP2, JPEG, JPG, PBM, PGM, PNG, PPM, TIF, TIFF, WBMP, WEBP, XBM, XPM and SVG image formats. You can export whole alignment or custom region. To select the custom region click on the *Select* button.