


Exporting Sequence Image

Use a sequence toolbar *Export image* button to save a screenshot of the sequence:



The Export Image dialog will appear where you should set the name, location, export settings, and format of the picture:

Export Image

Export settings

Area

- ☒ Currently viewed
- ☐ Zoomed annotations
- ☐ Sequence details

Export to file

File name

...

Format

PNG

▼

Help

Cancel

Export

UGENE supports export to the BMP, JPEG, JPG, PNG, PPM, TIF, TIFF, XBM, XPM, and SVG image formats. You can export currently viewed, zoomed annotations or sequence details areas. Also you can export whole sequence or custom region. Use the *Region* settings to do it.