


# Exporting Alignment as Image

To export an alignment as an image click the *Export as image* button on the editor toolbar or call the *Export->Export as image* context menu item.



The Export Image dialog will appear where you should set name, location, export settings and format of the picture:

A screenshot of the 'Export Image' dialog box. The dialog has a title bar 'Export Image'. It contains two main sections: 'Export settings' and 'Export to file'. In the 'Export settings' section, there is a 'Region' dropdown menu set to 'Whole alignment' with a 'Select' button next to it. Below this are three checkboxes: 'Include sequences' names' (unchecked), 'Include consensus' (unchecked), and 'Include ruler' (checked). In the 'Export to file' section, there is a 'File name' text field with a browse button ('...') to its right, and a 'Format' dropdown menu set to 'PNG'. At the bottom of the dialog are three buttons: 'Help', 'Cancel', and 'Export'.

UGENE supports export to the BMP, CUR, DDS, ICNS, JP2, JPEG, JPG, PBM, PGM, PNG, PPM, TIF, TIFF, WBMP, WEBP, XBM, XPM and SVG image formats. You can export whole alignment or custom region. To select the custom region click on the *Select* button.