


# Creating and Editing Qualifier

To add a qualifier to an annotation select it in one of the *Sequence View* subviews and press the Insert key, or use the *Add Qualifier* context menu or the *Annotations* main menu item.



The dialog will appear:

Add New Qualifier

Name

new\_qualifier

Value

Help

Cancel

OK

Here you can specify the name and the value of the qualifier.

To edit a qualifier, select the qualifier and press the F2 key or use the *Edit-->Qualifier* context menu item.

The dialog will appear:

**Edit Qualifier**

Name

Value