

Export Chromatogram Consensus

To export consensus sequence use the *Consensus* tab of the *Options Panel*:

The screenshot shows the 'Consensus' tab of the 'Options Panel'. The panel is divided into two main sections. The top section, 'Consensus mode', includes a 'Consensus type' dropdown menu set to 'Strict' and a 'Threshold' slider set to 100%. The bottom section, 'Export consensus', is highlighted with a red background and contains an 'Export to file' text box with the value 'Aligned reads_consensus_3.txt', a 'File format' dropdown menu set to 'Plain text', a 'Keep gaps' checkbox, and an 'Export' button.

The following parameters are available:

Export to file - here you need to select path for the output file.

File format - format for the output file.

Keep gaps - set this parameter if the result consensus must keep gaps.

When you click on the *Export* button the consensus sequence will be exported into selected output file.