


Exporting Sequence Image

Use a sequence toolbar *Export image* button to save a screenshot of the sequence:



The Export Image dialog will appear where you should set the name, location, export settings, and format of the picture:



You can export *Currently viewed*, *Zoomed annotations* or *Sequence details* areas. For zoomed annotations and sequence details areas you can export the whole sequence or custom region. Use the *Region* settings to do it.

Here you can browse for the file name and image format: JPG, BMP, SVG, PNG and TIFF.

Pay attention to *Quality* slider, it is present for JPG format only. Note, that SVG format is not supported for *Currently viewed area*.