


# Creating Fragment

To create a DNA fragment from a sequence region activate the *Sequence View* window and select either the *Actions Cloning Create Fragment* item in the main menu or the *Cloning Create Fragment* item in the context menu.

The *Create DNA Fragment* dialog appears:



If a region has been selected you can choose to create the fragment from this region. Otherwise you can either choose to create the fragment from the whole sequence or choose the *Custom* item and input the custom region.

To add a 5' overhang to the direct strand check the *Include Left Overhang* check box and input the required nucleotides. To add a 5' overhang to the reverse strand in addition to the described steps select the *Reverse-complement* item in the same group box.

Similarly, to add a 3' overhang check the *Include Right Overhang* check box, input the required overhang and select either the direct or the reverse-complement strand.

On the *Output* tab of the dialog you can optionally modify the annotations output settings.

Finally, press the *OK* button to create the fragment. The fragment will be saved as an annotation.