

Copying Qualifier Text

Use the Copy/Paste->Copy qualifier '[the qualifier name]' value context menu item to copy the qualifier value.

Auto-annotations [sars.gb | NC_004718]

NC_004718 features [sars.gb] *

> 3'UTR (0, 1)

> 5'UTR (0, 1)

▼ CDS (0, 14)

▼ CDS

codon_start

db_xref

db_xref

gene

locus_tag

note

product

protein_id

translation

> CDS

> CDS

> CDS

> CDS

> CDS

> CDS

> CDS

> CDS

> CDS

> CDS

> CDS

CDS

join(265..13398,13398..21485)

1

GI:30124074

GI:30124074

Find qualifier...

Invert annotation selection

Add 'db_xref' column

Go to position... Ctrl+G

Select sequence region... Ctrl+A

New annotation... Ctrl+N

Copy/Paste

Select

Add

Analyze

Align

Cloning

Export

Edit

Remove

Copy selected sequence Ctrl+C

Copy selected complementary 5'-3' sequence Ctrl+Shift+C

Copy amino acids Ctrl+T

Copy amino acids of complementary 5'-3' strand Ctrl+Shift+T

Copy annotation sequence

Copy annotation amino acids

Copy qualifier 'db_xref' value

Paste sequence Ctrl+V